



Silver Shadow

When eyewear becomes a fashion statement



Silver Shadow

Flash mirror on corrective lenses

HIGH PERFORMANCE WITH COSMETIC APPEAL... BY ESSILOR!

Silver Shadow by Essilor:

- A transparent mirror coating that does not alter the original tint ¹.
- A style that will not go unnoticed, offering a look that can be sporty, casual or elegant.
- An effect available on a range of tints and densities.



A HARD WEARING MIRROR COATING FOR A LONG-LASTING EFFECT

Silver Shadow is the product of Essilor's expertise in vacuum coating:



A LENS FOR EVERY TYPE OF CORRECTION... AND EVERY STYLE!

- Available across the entire range of single-vision and varifocal corrective lenses in Orma[®] and Styliis 1.67^{®3}.
- Anti-reflective coating on the concave side as standard to ensure optimal visual comfort.



¹ Silver Shadow does not alter the ISO sun-protection class of tinted lenses, with the exception of Class 0 lenses (tinted or not), which become Class 1 with Silver Shadow

² Performance equivalent to the Crizal[®] coating



















³ Silver Shadow is not available on photochromic and polarising lenses.

Silver Shadow can be combined with the most suitable tints on Orma® and Styli 1.67® single-vision and varifocal lenses

TINTS TRUE TO LIFE... PHYSIOTINTS®

Silver Shadow PHYSIOTINTS® means:







- A scientifically-developed range of tints offering wearers comfortable and natural colour.
- A natural perception of colour with a touch of originality.

	CLASS BEFORE SILVER SHADOW	CLASS WITH SILVER SHADOW	PHYSIO BROWN	PHYSIO GREY	PHYSIO GREY-GREEN	TRANSMISSION VALUE WITH SILVER SHADOW
Full tints	3	3				13%
Full tints	2	2				33%
Full tints	1	1				56%
Full tints	0	1				71%
Graduated tints	3	3				13%
Graduated tints	2	2				33%

FASHION TINTS... BEAUTY EYES®

- A range of delicate shades that bring out the wearers personality and gaze.
- Ideally suited to the image-conscious woman.

A combination of Silver Shadow and BEAUTY EYES® guarantee a stunning look for any wearer.

	LAVENDER	BLUEBELL	AMETHYST
Full tints 1			
Graduated tints 1			
MAXIMUM TRANSMISSION VALUE WITH SILVER SHADOW	65%	68%	60%

NATURALLY MYSTERIOUS... CLEAR LENSES

- Highly fashionable flash mirrors on non-tinted lenses play with the light giving the wearer a mysterious gaze.
- Turning everyday eyewear into a fashion statement.

	CLASS BEFORE SILVER SHADOW	CLASS WITH SILVER SHADOW		TRANSMISSION VALUE WITH SILVER SHADOW
Clear	0	0		82%

Silver Shadow

Technical performance that makes the difference!

Coatings included:

- Hard coats on both sides of the lenses
- Anti-reflection coating on the concave side
- Hydrophobic top coat on both sides
- Flash mirror on the convex side

Material characteristics:

Orma®

- Index (Ne): 1.502
- Density (g/cm³): 1.32
- Abbe number (√e): 59

Stylis®

- Index (Ne): 1.665
- Density (g/cm³): 1.36
- Abbe number (√e): 32

ISO TINT CLASSIFICATION AND USAGE RECOMMENDATIONS

	FILTER CLASS	CONDITIONS	USAGE
CLEAR	0		Indoors / Overcast
	1		Low sunlight Not suitable for driving at night
	2		Medium sunlight Not suitable for driving at night
	3		Bright sunlight Not suitable for driving at night
DARK	4		Exceptionally bright sunlight Not suitable for driving

----- Density not available with Silver Shadow. Information provided for reference purposes only -----

Issue Code 03.04

ESSILOR LIMITED, COOPER ROAD, THORNBURY, BRISTOL BS35 3UW
WWW.ESSILOR.CO.UK
 GENERAL ENQUIRIES 01454 417100

57-02-080

