



# Essilor Transitions® VI

## The New Range




**Transitions® VI**  
ADVANCED PERFORMANCE  
FOR HEALTHY SIGHT





## Essilor Transitions® VI

Offering the most advanced photochromic lenses with enhanced performance for healthy sight. Essilor Transitions® VI lenses automatically adapt to changing light conditions to reduce eye fatigue and eye strain and are darker outdoors particularly in hot temperatures.



**Transitions® VI gets darker in hot temperatures**


Transitions® VI lenses enhance the visual quality and help preserve the health and wellbeing of your patients' eyes, so they will see better today and tomorrow. Offering improved photochromic performance, Transitions® VI:

- > Is even darker outdoors, particularly in hot temperatures, for enhanced visual comfort in the sun
- > Is as clear as a clear lens indoors
- > Blocks 100% of UVA and UVB radiation
- > Performs consistently across all materials

**TRANSITIONS® VI IS 10% DARKER THAN TRANSITIONS® V IN HOT TEMPERATURES.**  
Darkness at 35°C/95°F Photochromic performance is influenced by temperature, UV exposure and lens material.



## Crizal® Alizé and Scotchgard™ enhances the performance of Essilor Transitions® VI which means the lenses return to clear faster



\* Transitions® VI with Crizal® Alizé and Scotchgard™ compared to Transitions® VI with hard coat.

Crizal® Alizé and Scotchgard™ enhances the performance of Essilor Transitions® lenses and offers the best resistance to smudges and dust. With the superior performance of Crizal® Alizé and Scotchgard™ your patients will appreciate:

- > Improved glare reduction, for enhanced visual comfort and performance, and better-looking lenses
- > Superior oleophobic and hydrophobic technology to minimise smudges, for a cleaner lens
- > The best antistatic properties to minimise dust attracted to the lens, so lenses remain cleaner for longer.

Available in all Single Vision and Bifocal as well as Varilux®

**TRANSITIONS® VI WITH CRIZAL® ALIZÉ AND SCOTCHGARD™ FADES BACK TO 70% TRANSMISSION 30% FASTER THAN TRANSITIONS® VI WITH ONLY A HARD COAT.**

# Essilor Transitions® VI lenses

are available in a complete range of designs and materials.

	Orma I.5	Airwear I.59	AS Airwear I.59	Ormix I.6	AS Ormix I.6	AS Stylis I.67	Available From
Single Vision Stock	+4.00/-4.00 cyls 2	+4.00/-6.00 cyls 2	n/a	+4.00/-8.00 cyls 2	n/a	n/a	1.10.08
Single Vision Rx	+7.00/-10.00 cyls 6	+6.00/-10.00 cyls 4	+8.00/-10.00 cyls 4	+8.00/-13.00 cyls 4	+10.00/-13.00 cyls 4	+8.00/-16.00 cyls 4	1.10.08
Single Vision 360°	n/a	n/a	+8.00/-10.00 cyls 4	+8.00/-13.00 cyls 4	n/a	n/a	1.01.09
Single Vision f-360°	n/a	n/a	+8.00/-10.00 cyls 4	+8.00/-13.00 cyls 4	n/a	n/a	1.01.09
Bifocal - 28 Flat Top	+5.00/-7.50 cyls 5	n/a	n/a	n/a	n/a	n/a	1.10.08

	Orma I.5	Airwear I.59	Ormix I.6	Stylis I.67	Available From
Varilux Ipeo	+6.00/-9.00 cyls 4	+6.00/-10.00 cyls 4	+6.00/-12.00 cyls 4	+8.00/-12.00 cyls 4	1.10.08
Varilux Physio	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 6	+8.00/-12.00 cyls 6	+9.00/-12.00 cyls 6	1.10.08
Varilux Physio f-360°	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 6	+8.00/-12.00 cyls 6	+9.00/-12.00 cyls 6	1.10.08
Varilux Panamic	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 4	n/a	+8.00/-12.00 cyls 4	1.12.08
Varilux Comfort	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 4	+8.00/-12.00 cyls 4	+8.00/-12.00 cyls 4	1.12.08
Varilux Comfort Advans	+6.00/-9.00 cyls 6	+6.00/-9.00 cyls 6	+8.00/-12.00 cyls 6	+8.00/-12.00 cyls 6	1.12.08
Varilux Ellipse	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 4	+8.00/-12.00 cyls 4	+8.00/-12.00 cyls 4	1.12.08
Varilux Ellipse Advans	+6.00/-10.00 cyls 6	+6.00/-10.00 cyls 6	+8.00/-12.00 cyls 6	+8.00/-12.00 cyls 6	1.12.08
Varilux Liberty	+6.00/-10.00 cyls 4	+6.00/-10.00 cyls 4	+6.00/-10.00 cyls 4	n/a	1.12.08

## COATINGS

Transitions® VI lenses come with a choice of coatings including Supra, Crizal® Alizé and the performance enhancing Crizal® Alizé and Scotchgard™.

## ORDERING

You can order Essilor Transitions® VI directly from Essilor or from a Varilux Specialist Laboratory or Wholesale Glazing Partner.

Essilor Limited  
Cooper Road  
Thornbury  
Bristol BS35 3UW  
General Enquiries: 01454 417100

[www.essilor.co.uk](http://www.essilor.co.uk)